



# FibroTest-ActiTest in HBV

Estimating Fibrosis and Activity - Detecting Inactive Carriers  
Recognized by EASL-ALEH<sup>1</sup>, APASL<sup>2</sup> and WHO<sup>13</sup> Guidelines



## Diagnosis

FibroTest offers the same highly accurate evaluation of both advanced fibrosis and cirrhosis regardless of ethnicity, HBeAg status, viral load, gender or ALT transaminases level.<sup>3</sup>

FibroTest evaluates both fibrosis in the early stages (F0 to F3) as well as cirrhosis (F4)<sup>3</sup>. Cirrhosis can be fine-tuned to 3 levels:<sup>5,6</sup>

- **F4.1** (cut-off 0.74): cirrhosis without complications
- **F4.2** (cut-off 0.85): cirrhosis with oesophageal varices only
- **F4.3** (cut-off 0.95): cirrhosis with severe complications (primary liver cancer; bleeding or decompensation)

## Follow-up

FibroTest and ActiTest may easily be repeated for the assessment of disease progression.<sup>6,7</sup>

FibroTest is validated both for the initial diagnosis of fibrosis and for the monitoring of patients, whether treated or untreated.<sup>6,7</sup>

ActiTest is more accurate than ALT transaminases for the diagnosis of necroinflammatory activity.<sup>8</sup>

## Identifying HBV Inactive Carriers

FibroTest F0  
AND  
ActiTest A0  
AND  
Viral Load low/undetectable

Fig 1: Identifying low-risk patients (HBV inactive carriers)

## Inactive carriers of HBV

- Using FibroTest and ActiTest in combination offers a more simple definition of HBV inactive carriers (Figure 1) <sup>4,9</sup>
- FibroTest has a superior prognostic value compared to both HBV viral load and ALT transaminases<sup>4,9</sup>

## FibroTest: the best option?

**vs. Transient Elastography (TE):** FibroTest identifies significant fibrosis with superior accuracy to TE, while identifying cirrhosis with similar accuracy<sup>3</sup> and much superior applicability without failure (98% vs 82%). Repeated FibroTests, unlike TE, do not have activity-related nor operator-related variability<sup>10,11</sup>.

**vs. APRI and FIB-4:** FibroTest does not include AST nor ALT transaminases, avoiding the risk of confounding features of fibrosis and activity<sup>12</sup>. For both cirrhosis and fibrosis, FibroTest is superior.<sup>3</sup>

### References:

1. EASL-ALEH Guidelines. J Hepatol 2015
2. Shiha G et al. Hepatol Int 2009
3. Houot M et al. Aliment Pharmacol Ther 2015 in press
4. Poynard T et al. J Hepatol 2014
5. Park MS et al. Liv Int 2015
6. Poynard T et al. Antivir Ther 2010
7. Haseltinej EL et al. Viral Hepat 2015
8. Poynard T et al. Gastroenterol Clin Biol 2010
9. Ngo Y et al. PLoS One 2008
10. Poynard T et al. Clin Chem 2010
11. Castera et al. Hepatology 2010
12. Poynard T et al. Clin Res Hepatol Gastroenterol 2014
13. Guidelines for the screening, care and treatment of persons with chronic hepatitis B infection, WHO, 2015

Assays (done at local lab.): Alpha-2 macroglobulin, Haptoglobin, Apolipoprotein A1, Total bilirubin, GGT, ALT, age, sex - according to BioPredictive precautions of use (biopredictive.com)

Find all the scientific publications of BioPredictive non-invasive tests on the website:

library.biopredictive.com

**BioPredictive S.A.**  
218 Boulevard Saint-Germain  
75007 PARIS - FRANCE  
Tel : +33 1 84 79 23 90  
contact@biopredictive.com

